

**MEDITERRANEAN  
PROTOCOL**

THE LASER REVOLUTION  
HAS A NAME:

*Red Edition*





# MEDITERRANEAN PROTOCOL

## Our revolution is your success

2024 is a crucial year for the future of professional aesthetics: It's the year you can enhance your practice with absolute confidence, thanks to the Mediterranean Protocol.

## The most revolutionary protocol to always be in the spotlight

Now you can work with effective fluence rates on skin types IV, V, and VI and on all hair types\*, from very light to the very dark, while maintaining safety throughout each treatment.

(\*except grey, white or red hair)

The innovative blend of wavelengths, scientifically proven to be effective and safe, reduces the risks of hyperpigmentation and burns. Additionally, our exclusive add-ons, such as the Asclepion Laser Skin Care X3<sup>®</sup>, significantly shorten treatment times. This advanced approach ensures high performance and safety for users.

*Red Edition*



# AN INNOVATIVE MIX

In our high-tech laboratories, we have developed a synergistic system of wavelengths to elevate laser hair removal treatments to a new level.

**The simultaneous emission of 808nm and 1060nm wavelengths** provides a wide range of absorption and penetration, allowing for a safe and effective treatment.



# POWERLINE MIX

## Where power meets sensuality

The Red Edition means combining the benefits of the Mediterranean Protocol with the advanced and leading features of a technology that has been leading the market for years.


It only takes a quick look at the basic technical data to highlight the performance of the PowerLine MIX handpiece, the international standard for performance and quality.

**POWER**  
**5,000 W**

**SPOT  
SIZE**

**3** cm<sup>2</sup>

**808** nm



**1,060** nm  
**WAVELENGTHS  
IN SYNERGY**

**SPEED** **20** hz

**PULSE DURATION**

**FROM** **3** ms



# THE KIT OF THE FUTURE FOR YOU

Only with Asclepion –  
Treatments without gel.  
Your starter kit to fully kickstart with  
the Mediterranean Protocol.



# NO MORE GEL!

Yes, that's right: The gel disappears forever from the treatments. We have replaced it with a specially developed product: **Asclepion Laser Skin Care X3®**.

## 1 CLEANSSES

By spraying it on the skin before starting the laser treatment, it thoroughly cleanses it

## 2 FACILITATES THE TREATMENT

It ensures smooth gliding of the handpiece without leaving sticky residues

## 3 REGENERATES

Excess product is gently massaged into the skin, providing excellent post-laser treatment care



**CONTAINS ALOE VERA**



## HYGIENE CAP

**Maximum comfort, maximum hygiene**

The disposable hygiene caps are the additional plus in customer care: Easy to use, they maximize hygiene and ensure clean treatments.

With the hygiene caps, you extend the lifespan of your handpiece!







# WELCOME TO THE RED EDITION BY ASCLEPION

**A pioneering marketing concept bringing together Asclepion's latest advancements in laser hair removal technology and treatment protocols for 2024.**

Under the Red Edition banner, all these breakthroughs are united by a consistent marketing approach.

Each product of the Red Edition has been developed and tested by a dedicated team of experts to meet the highest quality standards.

Joining the Red Edition means embracing cutting-edge technology and top-quality standards in the beauty industry.





# THE BENEFITS OF THE RED EDITION



All accessory packaging is designed to the smallest detail and tailored to the image of the **Red Edition.**



Our **dedicated marketing tools leverage** the Mediterranean Protocol to captivate your customers.

# WE ARE READY. EXPERIENCE WITH US A NEW ERA OF LASER HAIR REMOVAL!


- 1** High safety for **darker skin types**:  
Perfect for treatments in summer.
- 2** **Treatments without gel**: Annoying residues are replaced by **refreshing care**.

**ASCLEPION IS CONSTANTLY BY YOUR SIDE TO SUPPORT YOU WITH A PLAN OF INITIATIVES FOR YOUR HAIR REMOVAL BUSINESS.**

- 3** Treats **thick and deep hair** with **unprecedented effectiveness**.
- 4** **Reduces treatment times**:  
Faster treatments and **extra comfort**.



*Red Edition*





**INNOVATION 2024**  
FOR LASER HAIR REMOVAL



*I can*  
**SENSE**  
*it*



[www.asclepion.com](http://www.asclepion.com)

 **MADE IN GERMANY**

THIS BROCHURE IS NOT INTENDED FOR THE U.S. MARKET

Asclepion Laser Technologies GmbH | Brüsselerstrasse 10 | 07747 | Jena | Germany